







INSTALLATION OPTIONS
a) Concrete screw anchor
b) drop-in anchor


Use plywood template 30013914 to ensure an accurate installation.



Main direction towards patient chair headrest.

-Due to limitations in printing and copying procedures this template may not accurately correspond to the actual product.
-All dimensions are approximate.
-All dimensions are subject to change without a prior notification.
-Planmecca Oy refers the right to change this template and installation scheme without a prior notification.



Description of changes		Sheet #		Format	Status
Name		1/1		A0	
Material / Date					
Finish					
General Tolerance					
PLANMECCA		Dimensions generally according to drawing, (except if not supplied)		Scale	Revision
Assemblage: 6		1:1		1:1	
P.O. 20050 Helsinki, Finland		Drawing number			
Tel. +358 20 7795 500		30016615 - A			
www.planmecca.com					
Changing part geometry without written permission from Planmecca Oy is prohibited.					
© Copyright Planmecca Oy - all rights reserved.					